

BAILEY RIB APPROXIMATOR 18 CM

Ref. IR03100



27/08/2024 05:34:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

BAILEY RIB APPROXIMATOR 18 CM



Ref. IR03100

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040029016	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	3.6
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	20.5
Packed unit / Weight	0.162	Bare Unit / Weight (kg)	0.162
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 05:34:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr