

## PERCY AMPUTATION RETRACTOR (DIAM 17 CM) REMOVABLE HANDLES



Ref. IR05100



27/08/2024 05:35:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

# PERCY AMPUTATION RETRACTOR (DIAM 17 CM) REMOVABLE HANDLES



Ref. IR05100

## Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040029030	HS Code	90189090
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	19	Bare Unit / Width (cm)	19
Packed unit / Height (cm)	2	Bare Unit / Height (cm)	2
Packed unit / Length (cm)	19	Bare Unit / Length (cm)	19
Packed unit / Weight	0.6	Bare Unit / Weight (kg)	0.599
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 27/08/2024 05:35:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)