

AMALGAM PLUGGERS, SERRATED END, BLACK, NO.3



Ref. IDFABLA3



 27/08/2024 00:18:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

AMALGAM PLUGGERS, SERRATED END, BLACK, NO.3



Ref. **IDFABLA3**

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040015569	HS Code	9018.9090
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	1	Bare Unit / Width (cm)	1.8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.7
Packed unit / Length (cm)	18	Bare Unit / Length (cm)	18
Packed unit / Weight	0.027	Bare Unit / Weight (kg)	0.027
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 00:18:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr