

SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE PEDIATRIC FINGER CLIP FOR MP3



Ref. 222261



30/08/2024 16:49:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE PEDIATRIC FINGER CLIP FOR MP3



Ref. **222261**

Description

NELCORR OXIMAX compatible spo2 sensor, pediatric for oximeter MASTER PALM 3

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030148	HS Code	90181990
Warranty (months)	12	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.094
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 16:49:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr