

SILICONE TUBE ROLLER Ø6X10MMX50M NEW ASPIRET / ASKIR / EMIVAC



Ref. BCASPI144



 29/08/2024 17:36:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SILICONE TUBE ROLLER Ø6X10MMX50M NEW ASPIRET / ASKIR / EMIVAC



Ref. **BCASPI144**

Description

Roll silicone tube diam 6mm / ext 10mm x 50 m

Features

CE Marking	CE0123	Medical device	No
Bar code	0000000000895	HS Code	90330090
Warranty (months)	12	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	41	Packed unit / Height (cm)	13
Packed unit / Length (cm)	41	Packed unit / Weight	3.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 17:36:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr