

DERMAL CURETTES 2MM BOX OF 10

Ref. CPCURE01



29/08/2024 20:28:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

DERMAL CURETTES 2MM BOX OF 10



Ref. CPCURE01

Medical device used in light skin surgery procedures

Description

The very tapered sharp edge ring allows tissue harvesting and provides good visibility of the operating field.

The ergonomic handle allows a good grip of the instrument and ensures a safe and precise gesture.

- Equipped with a

- cap to prevent any risk of cut

- Sterilization by gamma irradiation 2.5 MEGA-RADS

- Packaging in individual sterile peelable bag covered with a transparent film

- Blue sleeve

- Length of the handle: 13 cm

- Box of 10

- Ready to use

- Complies with European directives 93/42/EC

Features

CE Marking	CE2460	Medical device	Class IIa
Bar code	3661040185996	HS Code	90183290
Warranty (months)	12	Sterile	Yes
Country	India		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	4	Packed unit / Height (cm)	7
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	13.2
Packed unit / Weight	0.13		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	10	Max temperature	40
Minimum humidity	35	Maximum humidity	65
OTHER FEATURES			
Packed unit / Quantity	1		

📅 29/08/2024 20:28:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr