


PUNCH BIOPSY Ø 4 MM BOX OF 10

Ref. CPPUNC07



29/08/2024 20:17:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PUNCH BIOPSY Ø 4 MM BOX OF 10



Ref. CPPUNC07

Ready-to-use, single-use device

Description

Efficient and precise thanks to its very sharp seamless stainless steel border. The notched neck ensures a good grip of the instrument. The size of the punch indicated on the handle allows easy identification. - Dimensions of the

handle: L 82 x diameter 10 mm

- Sterilization by gamma irradiation higher than the standard of 2.5 MEGA-RADS

- Under individual sterile peelable packaging covered with a transparent film

- Box of 10

- Complies with European directives 93/42/EC

Features

CE Marking	CE2460	Medical device	Class IIa
Bar code	3661040186214	HS Code	90183290
Warranty (months)	12	Sterile	Yes
Country	India		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	8.2
Packed unit / Weight	0.07		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	10	Max temperature	40
Minimum humidity	35	Maximum humidity	65

OTHER FEATURES

📅 29/08/2024 20:17:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr