

LEATHER DOCTOR'S ELITE BAGS BROWN



Ref. MEDOCT04_M



29/08/2024 17:48:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LEATHER DOCTOR'S ELITE BAGS BROWN



Ref. MEDOCT04_M

Essential and professional element of the work of the general practitioner in the context of his home visits and emergencies

Description

- Key closure- 3
- extractable compartments with transparent front part
- Isothermal amp maker, 31 elastics
- 1 transparent padded extractable compartment with Velcro separators for fragile instruments
- Mini-fastener
- Briefcase
- Shoulder strap and handle
- Material: leather
- Weight: 3.35 kg
- Outdoor measures: 41 x 31 x 14 cm / 17.8 L
- Recommended max. capacity: 4 kg
- Optional trolley
- Sold empty

Features

CE Marking	NONE	Medical device	No
Bar code	8436564540652	HS Code	4202.11.10
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	31
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	14
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	41
Packed unit / Weight	3.35	Bare Unit / Weight (kg)	4
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES: [3121-1-146-43](#) - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr