

INSTRUMENT TROLLEY S.S. 2 SHELVES WITH SIDE RAILS 60X40X80 CM



Ref. OSA6040



📅 09/10/2024 21:57:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

INSTRUMENT TROLLEY S.S. 2 SHELVES WITH SIDE RAILS 60X40X80 CM



Ref. OSA6040

High quality pedestal table, perfect for storing the necessary surgical instruments or consumables

Description

- Round tube structure in non-magnetic stainless steel for a better grip
- Stainless steel trays
- Removable galleries- Silent rubber wheels
- Permissible load: 80 kg
- Wheel diameter: 75 mm with brakes

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040168289	HS Code	9402900000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	45	Bare Unit / Width (cm)	40
Packed unit / Height (cm)	11	Bare Unit / Height (cm)	80
Packed unit / Length (cm)	85	Bare Unit / Length (cm)	60
Packed unit / Weight	8.6	Bare Unit / Weight (kg)	8.6
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 09/10/2024 21:57:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr