

LIAN NANO HAND-HELD MANOMETER/BLACK CASE/SILVER RING



Ref. 310065



CE
0459



30/08/2024 17:21:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LIAN NANO HAND-HELD MANOMETER/BLACK CASE/SILVER RING



Ref. 310065

Description

Dials 0 to 300 mmHg

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446023805	HS Code	90189010
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.423
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 17:21:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr