

## ASPIRATION PROBE CH 20 (10 PCS)



Ref. BCASPI200



 29/08/2024 16:17:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## ASPIRATION PROBE CH 20 (10 PCS)



Ref. BCASPI200

### Description

Tracheobronchial suction probes for CA-MI phlegm aspirators  
 Probes for use with CA-MI electric Suction Units: Aspiret / Askir 20 / Askir 30 / Askir C30 / Askir 36 BR / New Askir 36 BR Li-ON / New Askir 230/12v BR / New Aspiret

Single Use Probes  
 Box of 10 pieces.

### Features

CE Marking	NONE	Medical device	No
Bar code	3661040162072	HS Code	90189084
Warranty (months)	12	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.2
Bare Unit / Weight (kg)	0.2		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	1	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 29/08/2024 16:17:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
 13290 Aix-en-Provence | France

🌐 [www.gsh-med.fr](http://www.gsh-med.fr)

☎ +33(0)4 42 90 31 31

✉ [info@gsh-med.fr](mailto:info@gsh-med.fr)