

## SILICONE TUBE ROLL OF Ø6X10 MM - 1 METER



Ref. BCASPI201



29/08/2024 16:17:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## SILICONE TUBE ROLL OF Ø6X10 MM - 1 METER



Ref. BCASPI201

### Description

1 roll of tube silicone diam 6 / ext 10mm x 1 m

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040166353	HS Code	4202121900
Warranty (months)	24	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1
Packed unit / Height (cm)	nc	Packed unit / Length (cm)	nc
Bare Unit / Length (cm)	100	Packed unit / Weight	0.2
Bare Unit / Weight (kg)	0.2		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

29/08/2024 16:17:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)