

## SILICONE TUBE ROLL OF Ø8X14 MM - 1 METER



Ref. BCASPI202



29/08/2024 16:18:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## SILICONE TUBE ROLL OF Ø8X14 MM - 1 METER



Ref. BCASPI202

### Features

|                                    |               |                           |          |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking                         | CE0123        | Medical device            | No       |
| Bar code                           | 3661040169132 | HS Code                   | 90330090 |
| Warranty (months)                  | 12            | Country                   | Italy    |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |          |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | 0.2      |
| Bare Unit / Weight (kg)            | 0.2           |                           |          |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |          |
| Minimum temperature                | -20           | Max temperature           | 45       |
| Minimum humidity                   | NA            | Maximum humidity          | 85       |
| <b>OTHER FEATURES</b>              |               |                           |          |
| Packed unit / Quantity             | 1             |                           |          |

29/08/2024 16:18:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)