

## ELITE BAGS ISOTHERM BAG/CASE - COOL



Ref. MECOOL03\_B



07/11/2024 20:25:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## ELITE BAGS ISOTHERM BAG/CASE - COOL



Ref. MECOOL03\_B

Bag for transporting samples under optimal conditions

### Description

- Main isothermal compartment
- Temperature indicator- 4 removable isothermal subcompartments, including: -
- 1 compartment for urine samples, with removable separators, and a capacity of 24 tubes of 42 x 26.5 x 9.5 cm - 3 compartments for tube racks
- Capacity 300 tubes
- Thermometer included
- These compartments are removable, have a handle, an identifier and a sponge to absorb flows.
- Transparent upper pocket for documents
- Transparent front pocket with closure
- Material: tarpaulin- Weight: 5.30 kg
- Outdoor measures: 44 x 29 x 39 cm / 49 L
- Recommended max. capacity: 10 kg
- Bag in accordance with UN 3373, authorizing the transport of biological substances class B
- Sold empty

### Features

CE Marking	NONE	Medical device	No
Bar code	8436564540355	HS Code	4202.12.99.90
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	30	Bare Unit / Width (cm)	29
Packed unit / Height (cm)	38	Bare Unit / Height (cm)	39
Packed unit / Length (cm)	44	Bare Unit / Length (cm)	44
Packed unit / Weight (kg)	5.1	Bare Unit / Weight (kg)	5.3

#### SPECIFIC STORAGE CONDITIONS

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)