

## ELITE BAGS DOCTOR, DOCTOR CASE - BLACK



Ref. MEDOCT03\_N



29/08/2024 16:31:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## ELITE BAGS DOCTOR, DOCTOR CASE - BLACK



Ref. MEDOCT03\_N

Essential and professional element of the work of the general practitioner, in the context of his home visits or emergencies

### Description

- Key closure- 3
- extractable compartments
- Isothermal amp maker with 31 elastics
- Extractable transparent padded compartment with Velcro dividers for fragile instruments
- Mini-fastener
- 
- Briefcase- Shoulder strap and handle
- Lanyard for positioning on a trolley- Material:
- 300D polyester high tenacity fabric
- Weight: 3.35 kg
- Outdoor measures: 40 x 30 x 14 cm / 20.4 L
- Recommended max. capacity: 4 kg
- Sold empty

### Features

CE Marking	NONE	Medical device	No
Bar code	8436564540010	HS Code	420212
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	40	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	16	Bare Unit / Height (cm)	14
Packed unit / Length (cm)	33	Bare Unit / Length (cm)	40
Packed unit / Weight	3.35	Bare Unit / Weight (kg)	3.35
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

09 83 31 31 33 - For information purposes only, can be subject to change

### OTHER FEATURES

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 22 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)