

PAIR OF SCREW TYPE EARTIPS/BLACK/FOR "TRIANGLE" STETHOSCOPE




Ref. 506920



 30/08/2024 18:15:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

PAIR OF SCREW TYPE EARTIPS/BLACK/FOR "TRIANGLE" STETHOSCOPE



Ref. **506920**

Description

Spare part (x2) for Cardio Prestige® Cardio Sonic® Pediasonic® Magister® Laubry® Pulse® & Dual Pulse® stethoscopes

Features

CE Marking	NONE	Medical device	No
Bar code	3700446009410	HS Code	90330000
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.005
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:15:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr