

## FLOATING DIAPHRAGM (SMALL) FOR CARDIOPRESTIGE II SPENGLER



Ref. 507977



 30/08/2024 18:55:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## FLOATING DIAPHRAGM (SMALL) FOR CARDIOPRESTIGE II SPENGLER



Ref. 507977

### Description

Spare part for SPENGLER Cardio Prestige® II stethoscopes.  
Multifrequency fiberglass membrane for pediatric pavilion.

### Features

CE Marking	NONE	Medical device	No
Bar code	3700446030360	HS Code	90181990
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.001
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 18:55:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)