

## LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)



Ref. 513229



CE  
0459

DISPOSITIF MEDICAL  
CLASSE  
**Im**  
MEDICAL DEVICE

GARANTIE  
**2**  
2 ANS  
2 YEARS  
GUARANTEE

28/11/2024 22:33:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)



Ref. 513229

Light and robust, the blood pressure monitor that combines style and performance.

### Description

MANOMETRE EN CONTROLE QUALITE : VERIFIEZ QUE VOTRE LOT EST ACCEPTE, AVERTIR POUR CONTROLE PAR ECHANTILLONNAGE

### Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446007195	HS Code	90189010
Warranty (months)	24	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	7
Packed unit / Length (cm)	19	Packed unit / Weight	0.526
Bare Unit / Weight (kg)	0.5		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

28/11/2024 22:33:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)