

## REPLACEMENT BULB FOR MOBI




Ref. 513719



 30/08/2024 19:12:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## REPLACEMENT BULB FOR MOBI



Ref. 513719

### Description

Spare part: Pear for MOBI blood pressure monitor from SPENGLER.  
The design of the SPENGLER pear follows the geometry of all hands.  
Its volume being higher than other models, the inflation power is increased.

### Features

CE Marking	NONE	Medical device	No
Bar code	3700446029272	HS Code	90189010
Warranty (months)	12	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.043
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 19:12:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)