

## MAXI +3 WALL MOUNTED WITH SET OF 3 CUFFS (S, M, L)



Ref. 522707



CE  
0123

DISPOSITIF MEDICAL  
CLASSE  
**Im**  
MEDICAL DEVICE

GARANTIE  
2 ANS  
**2**  
2 YEARS  
GUARANTEE

30/08/2024 19:39:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## MAXI +3 WALL MOUNTED WITH SET OF 3 CUFFS (S, M, L)



Ref. 522707

**SPENGLER technology in a blood pressure monitor with a dial with a diameter of 147mm.**

### Description

HIGH MEASUREMENT  
 ACCURACY• Double tubing for linear pressure rise

COMFORT• Ambidextrous•  
 Remotely readable measurement•  
 2.5-meter spiral tubing• Large  
 volume ergonomic pear for quick inflation•  
 Storage basket behind the dial•  
 Adjustable wall mount

ARMBAND• Supplied with kit of 3 cuffs (S, M, L), nylon

### Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3700446031152	HS Code	90189010
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.154
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 19:39:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
 13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)