

GREY SPENGLER (XL) CUFF (32-48CM) FOR TENSONIC & AUTOTENSIO



Ref. 527465


CE



30/08/2024 19:46:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

GREY SPENGLER (XL) CUFF (32-48CM) FOR TENSONIC & AUTOTENSIO



Ref. 527465

Description

Extra large grey armband adult (xl) (32-48cm) for TENSONIC and AUTOTENSIO electronic blood pressure monitor from SPENGLER

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446041915	HS Code	90262040
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	3
Packed unit / Length (cm)	16	Packed unit / Weight	0.13
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 19:46:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr