

EXTRA LARGE ADULT CUFF 40-48CM FOR SPOT-CHECK 2VS



Ref. 546021



30/08/2024 19:53:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

EXTRA LARGE ADULT CUFF 40-48CM FOR SPOT-CHECK 2VS



Ref. 546021

Description

Extra large armbands (40-48cm) for SPENGLER 2VS spot-check station.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029685	HS Code	9018901000
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.102
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 19:53:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr