

CRYO PROFESSIONAL 170 ML - 60 APPLICATIONS 2 MM

Ref. DUCRYO02



29/08/2024 19:47:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CRYO PROFESSIONAL 170 ML - 60 APPLICATIONS 2 MM



Ref. DUCRY002

Effectively treats benign skin lesions

Description

- Cryotherapeutic agent, allowing to deliver a temperature of -55 ° C in less than 15 seconds
- Application time between 15 and 40 seconds
- Complete set including:
 - A 170 ml aerosol can
 - 2mm and / or 5 mm applicator sticks, usable depending on the different lesions
 - An instruction manual
- The sale of this medical equipment is reserved for health institutions, doctors and doctors' offices

Features

CE Marking	CE0344	Medical device	Class IIa
Bar code	8717484005811	HS Code	90189084
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	14
Packed unit / Length (cm)	14	Packed unit / Weight	0.3
SPECIFIC STORAGE CONDITIONS			
Max temperature	25		
OTHER FEATURES			

29/08/2024 19:47:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr