

LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)



Ref. 513210



CE
0459


DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

 29/11/2024 00:14:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)



Ref. 513210

Light and robust, the blood pressure monitor that combines style and performance.

Description

MANOMETRE EN CONTROLE QUALITE : VERIFIEZ QUE VOTRE LOT EST ACCEPTE, AVERTIR POUR CONTROLE PAR ECHANTILLONNAGE

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446009953	HS Code	90189010
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	7
Packed unit / Length (cm)	19	Packed unit / Weight	0.5
Bare Unit / Weight (kg)	0.4		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/11/2024 00:14:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr