

## SUNTECH TANGO® M2 MONITOR WITHOUT ECG (CONNECTED TO THE STRESS TEST)

# SUNTECH

Ref. DSU230700

CE  
0413

DISPOSITIF MEDICAL  
CLASSE  
**IIa**  
MEDICAL DEVICE

GARANTIE  
**2**  
ANS  
2 YEARS  
GUARANTEE



30/08/2024 22:46:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

# SUNTECH TANGO® M2 MONITOR WITHOUT ECG (CONNECTED TO THE STRESS TEST)



Ref. DSU230700

**Automatic blood pressure monitor for cardiac stress tests**

## Description

- Reliable automated blood pressure (BP) to reduce the risk of transcription errors and focus on the patient
- The Tango® M2 automatic PNI monitor features DKA technology developed specifically for the cardiac stress testing environment. With advanced proprietary algorithms, SunTech Medical is the only company in the market for motion-tolerant BP measurement.
- 2 Blood pressure measurement modes: with or without effort
- Allows multiple repeated automated BP measurements at short intervals.
- Use the auscultatory method with a microphone in the cuff to measure voltage and see Korotkoff noises on the screen.
- 7-inch LCD screen
- USB port for on-site upgrades

The device needs an ECG signal to optimize the measurement of this voltage. The Tango M2 is therefore offered with its autonomous ECG system or connected directly to the ECG of the stress test (GE, Mortara...).


## Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	10840935105701	HS Code	9018.90.5040
Warranty (months)	24		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	17.4
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	11.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	24
Packed unit / Weight	nc	Bare Unit / Weight (kg)	1.68
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 30/08/2024 22:46:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)