

SET OF 4 PAPER ROLLS FOR SPENGLER 4VS MONITOR PRINTER



Ref. 551643



 30/08/2024 20:06:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SET OF 4 PAPER ROLLS FOR SPENGLER 4VS MONITOR PRINTER



Ref. 551643

Description

Set of 4 rolls of SPENGLER 4VS Vital Signs Monitor Printer Paper

Features

CE Marking	NONE	Medical device	No
Bar code	3700446033125	HS Code	85049099
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.228
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:06:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr