

## LED FIBER OPTIC OTOSCOPE - SMARTLED 5500 PEDIATRIC



Ref. 570517P



30/08/2024 20:13:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## LED FIBER OPTIC OTOSCOPE - SMARTLED 5500 PEDIATRIC



Ref. 570517P

The SMARTLED 5500 otoscope has been developed with a lens to allow exceptional diagnosis.

### Description

- Cutting-edge LED technology.
- Powerful luminance of 50K lumens,
- 4000K colorimetry, the perfect range for observation,
- CRI (Color Rendering Index) of 95 out of 100 for perfect color restitution.
- Optical fiber for direct, homogeneous and glare-free illumination of the ear canal and eardrum

Magnifying glass x3

#### EASY TO USE AND PRACTICAL

- Ergonomic handle
- ABS and stainless steel head for secure fixation
- Fixing clip with switch
- Crimped bulb for greater reliability
- LED bulb with a lifespan of 55000 h

#### ADDITIONAL INFORMATION

Comes with rigid housing consisting of 2 LR6 AA batteries, 8 speculums (2.5mm and 4mm)

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446041694	HS Code	90189042
Warranty (months)	36	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	12	Packed unit / Height (cm)	4
Packed unit / Length (cm)	20	Packed unit / Weight	0.385
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	90
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 20:13:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)