

## GUEDEL CANNULA SINGLE USE STERILE ISO 6.5 (WHITE/N°1)



Ref. 580854

CE  
0120

DISPOSITIF MEDICAL  
CLASSE  
**I**  
MEDICAL DEVICE

STERILE



16/09/2024 10:12:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## GUEDEL CANNULA SINGLE USE STERILE ISO 6.5 (WHITE/N°1)



Ref. 580854

### Description

Single-use cannula in individual sterile bag.  
The Guedel cannula is used to keep the airway of an unconscious patient open.

### Features

CE Marking	CE0120	Medical device	Class I
Bar code	3661040215891	HS Code	90183900
Sterile	Yes	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.008
Bare Unit / Weight (kg)	0.008		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

16/09/2024 10:12:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)