

LANCET 0,32MM 30G (200 UNITS)

Ref. 745013



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2
2 YEARS
GUARANTEE

STERILE

30/08/2024 20:30:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LANCET 0,32MM 30G (200 UNITS)



Ref. 745013

Description

Lances for use with the TB 100 Lancing Pen.

- Needle diameter 30G = 0.298mm ~ 0.320mm
 -Sterile, single use - Length 32mm +/- 1mm- Width 6.5mm +/-

0.5mm
 -Outer needle 3.3mm +/- 0.3mm- Packaging: box of 200 (

20 sachets of 10 lancets)

Find our complete kit Glucometer SPENGLER TB 100 (ref. 745000)

Refund LPP code: 6111373

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446041687	HS Code	90183190
Warranty (months)	24	Sterile	Yes
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	5
Packed unit / Length (cm)	10	Packed unit / Weight	0.22
SPECIFIC STORAGE CONDITIONS			
OTHER FEATURES			
Code LPP	6111373	Packed unit / Quantity	1

30/08/2024 20:30:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr