

## SPO2 EAR CLIP SENSOR FOR 4VS MONITOR



Ref. 551636


CE  
0537



 30/08/2024 20:05:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## SPO2 EAR CLIP SENSOR FOR 4VS MONITOR



Ref. 551636

### Description

SPO2 ear sensor for SPENGLER 4 VS vital signs monitor

### Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446031244	HS Code	90181990
Country	Uruguay		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.082
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 20:05:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)