

RECHARGEABLE BATTERY FOR SPENGLER 4VS MONITOR



Ref. 551641

 30/08/2024 20:05:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

RECHARGEABLE BATTERY FOR SPENGLER 4VS MONITOR



Ref. 551641

Description

SPENGLER 4VS Vital Signs Monitor Rechargeable Battery

Features

CE Marking	NONE	Medical device	No
Bar code	3700446033101	HS Code	85049099
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.147
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:05:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr