

SET OF SILICONE TUBES FOR AUTOCLAVABLE JARS FOR NEW ASPIRET, ASKIR



Ref. BCASPI136



29/08/2024 17:35:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SET OF SILICONE TUBES FOR AUTOCLAVABLE JARS FOR NEW ASPIRET, ASKIR



Ref. **BCASPI136**

Description

Set of pipes containing:- 1 pipe diam 6 mm / ext 10mm x length

140cm- 1

pipe diam 6 mm /ext 10mm x length 13 cm

- Connec connectors 8-9-10 mm

- 1 bacterial filter

Features

CE Marking	CE0123	Medical device	No
Bar code	3661040163604	HS Code	90330090
Warranty (months)	12	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	4.34
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 17:35:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr