

CA-MI NEW ASKIR 30 PROXIMITY-JAR FLOVAC 1L SUCTION UNIT



Ref. BCFLOV404



29/08/2024 21:20:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CA-MI NEW ASKIR 30 PROXIMITY-JAR FLOVAC 1L SUCTION UNIT



Ref. BCFLOV404

Fixed electric surgical Suction Unit, intended to drain body fluids and to perform oral, nasal and tracheal aspiration in adults and children.

Description

Suitable for tracheostomy patients, ideal for minor surgical applications as well as postoperative home treatments. The "PROXIMITY" function is used to activate or deactivate the device without touch contacts and prevents cross-contamination between two patients.

Main applications: Endoscopy, Gastroenterology, General surgery, Home care, Obstetrics, Room, Veterinary Med.

Composed of: New Askir 30 Proximity Suction Unit + Flovac 1L Kit (1x Flovac 1L reusable jar / 1x Sterile UU pocket 1L Flovac / 1x Set tubes Ø6X10mm + connectors / 1x Round Retractor for Flovac System / 1x Yellow Male Connector FLOVAC System (5 pcs) / 1x Configuration Manual Suction Units

Oil-free piston pump requiring no maintenance-

Main system is equipped with a pressure gauge (bar & kPa) and a vacuum regulator

- Supplied accessories: Fluid collection jar with an anti-reflux valve with different options / 6 x 10 Silicone hoses (autoclavable) and conical connectors ø 8-9-10 mm / Antibacterial & hydrophobic filter (single patient) / CH20 cannula / Power cord with Schuko plug

Power supply: 230V-50Hz - other configurations available on request

Maximum vacuum (adjustable): 0.80 Bar -80 kPa -600 mmHg-

Max free air flow: 40 l/min

Noise level: 60.5 dBA

Uninterrupted operating cycle


Autoclavable Makrolon polycarbonate jar

Flow regulation button and vacuum gauge

 29/08/2024 21:20:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

CA-MI NEW ASKIR 30 PROXIMITY-JAR FLOVAC 1L SUCTION UNIT



Ref. BCFLOV404

Features

CE Marking	CE0123	Medical device	Class IIa
HS Code	90189084	Warranty (months)	24
Country	Italy		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3.6
Bare Unit / Weight (kg)	3.6		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			

29/08/2024 21:20:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr