

## CA-MI NEW ASKIR 30 PROXIMITY-JAR FLOVAC 2L SUCTION UNIT



Ref. BCFLOV405



29/08/2024 21:20:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CA-MI NEW ASKIR 30 PROXIMITY-JAR FLOVAC 2L SUCTION UNIT



Ref. BCFLOV405

Fixed electric surgical Suction Unit, intended to drain body fluids and to perform oral, nasal and tracheal aspiration in adults and children.

### Description

CA-MI New Askir 30 Proximity Suction Unit - FLOVAC® 2L

Jar Fixed electric surgical Suction Unit, intended to drain body fluids and to practice oral, nasal and tracheal aspiration in adults and children.

Suitable for tracheostomy patients, ideal for minor surgical applications as well as postoperative home treatments. The "PROXIMITY" function is used to activate or deactivate the device without touch contacts and prevents cross-contamination between two patients.

Main applications: Endoscopy, Gastroenterology, General surgery, Home care, Obstetrics, Room, Veterinary Med.

Composed of: New Askir 30 proximity Suction Unit + Flovac 2L Kit (1x Reusable Flovac 2L jar / 1x Sterile UU pocket 2L Flovac / 1x Set tubes Ø6X10mm + connectors / 1x Round Retractor for Flovac System / 1x Yellow Male Connector FLOVAC System (5 pcs) / 1x Configuration Manual Suction Units

}-

Oil-free piston pump requiring no maintenance-

Main system is equipped with a pressure gauge (bar & kPa) and a vacuum regulator

- Supplied accessories: Fluid collection jar with an anti-reflux valve with different options / 6 x 10 Silicone hoses (autoclavable) and conical connectors ø 8-9-10 mm / Antibacterial & hydrophobic filter (single patient) / CH20 cannula / Power cord with Schuko plug

Power supply: 230V-50Hz - other configurations available on request

Maximum vacuum (adjustable): 0.80 Bar -80 kPa -600 mmHg-

Max free air flow: 40 l/min

Noise level: 60.5 dBA

Uninterrupted operating cycle

 29/08/2024 21:20:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

Antibacterial filter

Dimensions: 300 x 21 x 100 mm / Volume: 3ka60

3/3

## CA-MI NEW ASKIR 30 PROXIMITY-JAR FLOVAC 2L SUCTION UNIT



Ref. BCFLOV405

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040200088	HS Code	9018908400
Warranty (months)	24	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3.6
Bare Unit / Weight (kg)	3.6		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			

29/08/2024 21:20:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)