

DIGITAL OTOSCOPE FIREFLY DE500A

Ref. DFIDE500A



30/08/2024 09:24:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

DIGITAL OTOSCOPE FIREFLY DE500A



Ref. DFIDE500A

Connected otoscope for telemedicine use

Description

High-precision digital handheld otoscope video; makes very precise and clear observations of the ear.

- Resolution of 1280x1024 real pixels without processing with interpolation techniques.
- Provides instant images and videos: DE500 uses a state-of-the-art image processing engine that virtually eliminates visual lag, allowing for an intuitive user experience.
- Magnification up to 50x optically and 150x digitally for accurate analysis.
- Equipped with multi-layered glass lenses, rather than plastic, thus offering superb depth of field and unparalleled image sharpness.
- High-quality ultra-bright LEDs with adjustable brightness are installed.
- Easily integrates with digital medical records (EMRs)

Features

CE Marking	CE	Medical device	Class I
Bar code	8512890070640	HS Code	90189042
Warranty (months)	12		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	20	Packed unit / Height (cm)	5
Packed unit / Length (cm)	15	Packed unit / Weight	0.302
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 09:24:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr