

## PROBE AMPULARIUM, RED



Ref. MEPROB02\_R



29/08/2024 16:43:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## PROBE AMPULARIUM, RED



Ref. **MEPROB02\_R**

### Description

Isothermal amp maker designed to offer the professional a large capacity tool for transporting bulbs.

This innovative, functional and handy ampoulier allows the professional to transport in an orderly way and in perfect conditions the drugs necessary to offer quality care.

- 3 totally adjustable elastics that facilitate the insertion of bulbs and various materials of any format. Capacity for 132 ampoules- 8 extractable mesh pockets for

the insertion of the cold accumulator, bulbs or utensils of usual use (4 extractable)

- 2 external flat pockets for the pocket container of bio-contaminated material and other utensils

- Retractable handles- Card holder- Material: Polyester 600 x

600

- Weight: 750 g

- Outdoor measures: 28 x 21 x 9 cm / 5.29 L

- Maximum recommended capacity: 3.5 kg

- Sold empty

### Features

CE Marking	NONE	Medical device	No
Bar code	8436564540119	HS Code	420212
Warranty (months)	24	Country	China

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	23	Bare Unit / Width (cm)	21
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	9
Packed unit / Length (cm)	31	Bare Unit / Length (cm)	28
Packed unit / Weight	1.122	Bare Unit / Weight (kg)	0.75

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

#### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

29/08/2024 16:43:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)