

EXAMINATION COUCH ELECTRIC CAIX WITH WHEELS, BLUE

Ref. OV004_B



29/08/2024 16:49:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

EXAMINATION COUCH ELECTRIC CAIX WITH WHEELS, BLUE



Ref. OV004_B

Electric table from 560 to 960 mm - Width 670mm - with casters

Description

- Length of sections: backrest 60cm, seat 135cm. Width 67cm.
- Maximum load: 220kg dynamic, 400kg static.
- Product delivered assembled, with a roller holder.
- Weight of the sofa 80 kg.
- Level of protection of electric motors class: CLASS 2.
- Waterproof level of the remote control: IP 66.

Structure (frame and base):

- Material and coating of the structure, steel quality 102, covered with white epoxy paint.
- Dimensions of the tube section: frame: 40 X 20 cm, thickness 2 mm, base: 60 X 30, thickness 2 mm and diam 50.
- Securing the MIG chassis with welded reinforcement brackets under the connecting bars (4 bars) located under the seat and at the front base and rear base.
- Adjustable feet.
- Optional retractable casters.
- Compatibility of the couch with all patient lifts.

Top :

- Material: beech, 10 ply, thickness 16mm.
- Mode of attachment of the tray on the frame: for the seat, 8 pulls of 10 x 60 mm; The backrest is connected to the seat by 2 hinges fixed with 8 inserts of 6 mm diameter.
- Electric variable height from 52 to 92 cm.
- Linear cylinder, stroke 200 mm, power 6000 W.
- Manual variable height control lever.
- Supply voltage and frequency: 230 V, 50/60 HZ.

Backrest :

- Inclination assisted by gas cylinder.
- Cable control with joystick attached to the frame to the right of the sofa.
- Maximum inclination + 75°.

Upholstery

- Upholstery: coated fabric, polyurethane coating on 100% cotton jersey support.
- PF polyurethane foam, thickness: 4 cm, density: 45 kg/m³.
- Gluing the pvc on the polyurethane foam then stapling under the seat.
- Backrest, seat and leggings in sewn edges, then stapled underneath.

Fire behaviour:

- PVC in accordance with standard NF EN 1021.1 and EN 1021.2, non-fire classification M1.
- Foams seated, backrest, leggings, headrest, classification non-fire M4.
- Wood (multi-ply beech), non-fire classification M4. 3-year warranty:

mechanical parts and engine

1-year warranty: coating

 29/08/2024 16:49:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

EXAMINATION COUCH ELECTRIC CAIX WITH WHEELS, BLUE



Ref. OV004_B

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040094403	HS Code	94029000
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	67
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	92
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	195
Packed unit / Weight	nc	Bare Unit / Weight (kg)	82.85
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 16:49:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr