

**MEDIPLUGS OPERATING ROOM  
SHOES, NON-PERFORATED, 35, NAVY  
BLUE**

**medi Plogs**

Ref. WSAM\_41L

CE

 16/09/2024 22:38:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE



Ref. WSAM\_41L

=====

### Description

Clogs suitable for modern hospitals, operating theatres and intensive care units.

Products designed and manufactured in collaboration with researchers and doctors.

- Antistatic clogs with non-slip and non-chafing soles
- The back of the heel is sloped and rounded to make walking easier
- Ventilation system allows air to circulate
- Toe prints to enhance hold and spurs to aid the circulation
- 100% CFC-free polyurethane
- Washable by hand or machine at 60°C
- Sizes: 35 - 47

### Features

Bar code	3661040100074	HS Code	64029993
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	22
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	11
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	31
Packed unit / Weight	1.16	Bare Unit / Weight (kg)	1
<b>SPECIFIC STORAGE CONDITIONS</b>			
<b>OTHER FEATURES</b>			
Colors	fuchsia	Packed unit / Quantity	1
Shoe Size	41		

16/09/2024 22:38:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)