

CARDIOMATE EVI 7" TABLET ECG SPENGLER



Ref. 260050

CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
3
3 ANS
GUARANTEE



13/11/2024 00:45:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CARDIOMATE EVI 7" TABLET ECG SPENGLER



Ref. 260050

Spengler expands its range of ECGs and offers you its latest innovation: the CARDIOMATE EVI 7" tablet, a fully autonomous technological product.

Description

Very easy to use:

- 7" color touch screen
- Portable: compact size and light weight, easy to carry: 0.8kg and 197 x 112.4 x 26.1 mm
- Large battery capacity, rechargeable on mains
- Allows to print ECGs directly on a printer via wifi network

High-performance technology:

- High accuracy: sampling rate 32,000 Hz
- Detection and marking of the pacemaker
- Automatic interpretation (Glasgow) and calculation of intervals
- Input of patient information using a keyboard or barcode reader
- real time visualization of the 12 tracks at the same time as automatic interpretation and interval calculations. (possibility of a 18 leads tablet optional)
- Data export via wifi or USB port under .pdf or .xml

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6971049341789	HS Code	9018110000
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	10	Packed unit / Height (cm)	28
Packed unit / Length (cm)	33	Packed unit / Weight	2.8
Bare Unit / Weight (kg)	0.8		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	95
OTHER FEATURES			

📅 13/11/2024 00:45:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

🌐 www.gsh-med.fr