

POWER CABLE FOR NEW ASKIR 30



Ref. BCASPI170



 29/08/2024 21:25:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

POWER CABLE FOR NEW ASKIR 30



Ref. BCASPI170

Features

CE Marking	CE0124	Medical device	No
Bar code	3661040031552	HS Code	4202121900
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.158
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	5	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 21:25:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr