

## FLOVAC SET 1L FOR ASKIR 20/30/PROXIMITY/230/ASPIRET



Ref. BCASPI500



29/08/2024 20:55:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## FLOVAC SET 1L FOR ASKIR 20/30/PROXIMITY/230/ASPIRET



Ref. BCASPI500

### Description

Flovac 1L kit, consisting of: 1x BCASPI102: Reusable Flovac 1L jar for Aspiret, Askir 20 and 30 - 31843 1x BCASPI99: Sterile UU pocket 1L Flovac Aspiret, Askir 20 30 230 12V,

36BR 31848 1x  
 BCASPI111: Set tubes Ø6X10mm + connectors Flovac ASKIR 20/30/230/ASPIRET 1x  
 BCASPI105: Round Retractor for Flovac System New Askir 20,30,230/12VBR 1x  
 BCASPI164: Yellow Male Connector FLOVAC System (5 pcs) sf Askir12V/Emivac  
 1x Configuration Manual Phlegm  
 Suction Units

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040191676	HS Code	4202121900
Warranty (months)	24	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

29/08/2024 20:55:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
 13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)