

THERMOSCAN PRO 6000 RECHARGEABLE BATTERY PACK 1

WelchAllyn®

Advancing Frontline Care™

Ref. DWPRO6K09



29/08/2024 20:59:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

THERMOSCAN PRO 6000 RECHARGEABLE BATTERY PACK 1

WelchAllyn®

Advancing Frontline Care™

Ref. DWPRO6K09

Description

- 2.4V 750mAh rechargeable battery for Pro 6000 thermometer
- Fully charged, the thermometer can make 700 measurements
- Must be recharged thanks to the charging dock


Features

CE Marking	CE	Medical device	No
Bar code	00732094325089	HS Code	850610
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	2	Packed unit / Height (cm)	1
Packed unit / Length (cm)	4	Packed unit / Weight	0.042
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 29/08/2024 20:59:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr