

ELECTRIC EXAMINATION COUCH ROYE BLACK 67 CM

Ref. OV002_N



 29/08/2024 21:46:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ELECTRIC EXAMINATION COUCH ROYE BLACK 67 CM



Ref. OV002_N

Electric table from 490 to 960 mm

Description

- Built in white rectangular epoxy tube.
- Protective housing at the base.
- Equipped with 4 adjustable black pads, polyamide.
- Stitched upholstery, thickness 40 mm.
- Adjustable backrest by gas cylinder (bilateral control).
- Roller holder.
- Foot control.
- Seat foam and backrest density 28 kg / m3.
- Class 2 engine.

Dimensions

- H. all: L. 2000 × L. 670 × H. 490/960 mm.

- Backrest: L. 600 × L. 670 mm.

- Seat: L. 1350 × L. 670 mm.

- Sofa weight: 80 kg.

- Supported weight: 220 kg.

- 3-year warranty:
mechanical parts and engine

- 1-year warranty: coating

 29/08/2024 21:46:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ELECTRIC EXAMINATION COUCH ROYE BLACK 67 CM



Ref. OV002_N

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040203201	HS Code	94029000
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	91	Bare Unit / Width (cm)	67
Packed unit / Height (cm)	61	Bare Unit / Height (cm)	49
Packed unit / Length (cm)	205	Bare Unit / Length (cm)	200
Packed unit / Weight	nc	Bare Unit / Weight (kg)	100
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 21:46:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr