



SEMI-AUTOMATIC DEFIBRILLATOR FRED PA-1 SCHILLER




Ref. SSDEFI04



 16/09/2024 22:59:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SEMI-AUTOMATIC DEFIBRILLATOR FRED PA-1 SCHILLER



Ref. SSDEFI04

Simple semi-automatic automated external defibrillator (AED) within everyone's reach

Description

The FRED® PA-1 is perfectly suited for use on both adult and child patients thanks to its dedicated electrodes.

Suitable for use in healthcare facilities, ERPs: it is by design a device intended for the general public.

- Orange button to deliver electric shock (for semi-automatic version)

- Simple and uncluttered user interface, speech synthesis in French

- Has various cardiac massage aid functions to guide the user to perform quality CPR at the right frequency- Instantly turned on as soon as the hood is opened- Equipped by default with a metronome allowing the user to know the adequate rhythm of cardiac massage

- Performs periodic self-tests on a daily or weekly basis

- Default in French (Trilingual option possible)

- Can record up to 2h of data and 500 events - Pair of pre-connected defibrillation electrodes housed in the AED cover with the first aid kit* (slot for a second pair of electrodes under the AED)

- Indicator lights visible at more than 15m - Medical grade ABS V0 (flame retardant) <

font style="vertical-align: inherit;">- Maximum operating sound level: 85 dB - Comes standard with 1 battery: LiMn02 15V, 2.8Ah

*First aid kit: 1 pair of emergency scissors, 1 disposable razor, 1 pair of vinyl gloves, 1 ventilation mask, 2 sterile compresses, 1 towel to wipe the victim

 16/09/2024 22:59:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SEMI-AUTOMATIC DEFIBRILLATOR FRED PA-1 SCHILLER



Ref. **SSDEFI04**

Features

CE Marking	CE0459	Medical device	Class IIb
Bar code	3661040218328	HS Code	90189084
Warranty (months)	120	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	40	Bare Unit / Width (cm)	11
Packed unit / Height (cm)	15	Bare Unit / Height (cm)	30
Packed unit / Length (cm)	40	Bare Unit / Length (cm)	25.5
Packed unit / Weight	nc	Bare Unit / Weight (kg)	3.155
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

16/09/2024 22:59:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr