

SEMI AUTOMATIC DEFIBRILLATOR FRED PA-1 SCHILLER



Réf. SSDEFI02

SCHILLER


The Art of Diagnostics



CE
0459



SEMI AUTOMATIC DEFIBRILLATOR FRED PA-1 SCHILLER



Réf. SSDEFI02

FEATURES			
Medical device	Yes	Classe IIb	
Marquage	CE0459	Barcode	3661040218328
Guarantee (month):	120	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	11,0	Width packaging (cm):	40,0
Length (cm) :	25,5	Length packaging (cm) :	40,0
Height (cm) :	30,0	Height packaging (cm) :	15,0
Weight (kg) :	3,155	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



01/08/2022 12:48 - This article is available on our website : www.gsh-med.fr with the following reference SSDEFI02