

## CHILD SENSOR "DEPORTED" VELCRO FOR OXIMETER PC-60E



Ref. DPOXY05



 25/08/2024 20:41:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## CHILD SENSOR "DEPORTED" VELCRO FOR OXIMETER PC-60E



Ref. DPOXY05

### Description

With its remote pediatric sensor, disposable VELCRO. It will adapt to the feet of infants from 0 to 2 years old. The PC 60E pulse oximeter is reliable and accurate. Ideal for urban healthcare professionals.

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040165998	HS Code	9018909911
Warranty (months)	6	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	15	Packed unit / Height (cm)	1
Packed unit / Length (cm)	14	Packed unit / Weight	0.059
Bare Unit / Weight (kg)	0.059		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

25/08/2024 20:41:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)