

MOHR TUBING CLAMP WITH SCREW 18X18



Ref. IS13618



27/08/2024 06:41:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MOHR TUBING CLAMP WITH SCREW 18X18



Ref. IS13618

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040030777	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1.8
Packed unit / Height (cm)	nc	Packed unit / Length (cm)	nc
Bare Unit / Length (cm)	1.8	Packed unit / Weight	0.015
Bare Unit / Weight (kg)	0.015		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 06:41:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr