

ELITE BAGS ISOTHERM BAG/CASE - ROW'S



Ref. MEROW03_B



📅 25/08/2024 18:37:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

ELITE BAGS ISOTHERM BAG/CASE - ROW'S



Ref. MEROW03_B

Portable insulated case of medium capacity for transporting biological samples.

Description

Designed for transport to hospitals or homes of biological samples and analysis tubes. Very practical thanks to its shoulder strap that allows it to be easily worn on the shoulder and its interior that protects and maintains the temperature of the samples. Medium size. - Weight: 1 kg- Dimensions: 26 x

- 26 x 16 cm
- Material: tarpaulin
- Color: blue
- Capacity: 10 L
- Recommended max. content: 5 kg
- Sold empty

Features

CE Marking	NONE	Medical device	No
Bar code	8436564543554	HS Code	4202129990
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Bare Unit / Width (cm)	16
Packed unit / Height (cm)	26	Bare Unit / Height (cm)	26
Packed unit / Length (cm)	26	Bare Unit / Length (cm)	26
Packed unit / Weight	nc	Bare Unit / Weight (kg)	1
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

25/08/2024 18:37:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr