

METAL FEMALE CATHETER 15 CM (SET OF 4 PCS DIAM 2,3,4,5 MM)



Ref. IS03405



27/08/2024 06:20:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

METAL FEMALE CATHETER 15 CM (SET OF 4 PCS DIAM 2,3,4,5 MM)



Ref. IS03405

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040087597	HS Code	9018908500
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.7
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16.3
Packed unit / Weight	0.011	Bare Unit / Weight (kg)	0.011
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 27/08/2024 06:20:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr