

TEMPO WATERPROOF DIGITAL THERMOMETER SPENGLER

Ref. THE001



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2
2 ANS
2 YEARS
GUARANTEE

17/09/2024 13:01:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

TEMPO WATERPROOF DIGITAL THERMOMETER SPENGLER

Ref. THE001

TEMPO is a simple and precise digital thermometer

Description

TEMPO is a waterproof digital thermometer allowing 3 different uses for taking temperature: the oral, the rectal and the ampit one.

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6970392217642	HS Code	9025199090
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.012		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	95
OTHER FEATURES			

17/09/2024 13:01:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr